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
## Long-read capture with Twist target enrichment system

### Long-read capture with Twist target enrichment system

Tina Han, Holly Corbitt, Leonardo Arbizu, Esteban Toro, Chad Locklear.

Twist Bioscience, South San Francisco, California, USA

#### 1. Introduction

Target enrichment allows for high-resolution characterization of gene regions at a scale not possible with whole genome sequencing. While long-read PacBio HiFi sequencing has been shown to accurately and comprehensively characterize complex genomic regions, these reads are often limited by low coverage and high error rates. We describe a method to leverage Twist Bioscience target enrichment technology for gene panels sequenced with HiFi reads. This method increases the number and the alignment rate of sequencing reads, resulting in higher coverage and lower error rates. The resulting panels are highly accurate and can be used for a wide range of applications.

#### 2. Methods

We designed gene panels of various target sizes, ranging from 0.2-22 Mb. Panels were optimized using a proprietary algorithm to create balanced panels of varying sizes and the same base frequency distribution of all target regions. Our long-read hybrid capture protocol starts with 200-1000x fragmented HiFi reads that were enriched using Twist Bioscience target enrichment technology. The enriched libraries were then sequenced on a HiFi sequencer. The resulting panels were then used to generate a panel of 49 clinically significant genes. A set of 49 genes was selected for target enrichment. The panel captures reads from target regions and non-target regions. The panel captures reads from target regions and non-target regions. The panel captures reads from target regions and non-target regions.

Twist Bioscience target enrichment technology was used to generate HiFi reads, remove PCR duplicates, and demultiplex, and a HiFi sequencing pipeline was used to call variants for individual samples. The Pacific Bioscience HiFi sequencing pipeline is publicly available on GitHub: <https://github.com/PacificBiosciences/HiFiSequencingPipeline>

#### 3. Long-Read Pharmacogenomics (PGx) Panel

Pharmacogenomics, including HLA genes and CYP2D6, are increasingly difficult to genotype via array or via short-read sequencing as the need to have accurate candidate and high-coverage sequencing. CYP2D6 is also difficult to genotype via array or via short-read sequencing as the need to have accurate candidate and high-coverage sequencing. CYP2D6 is also difficult to genotype via array or via short-read sequencing as the need to have accurate candidate and high-coverage sequencing.

Panel of 49 Clinically Significant Genes (full gene coverage in bold)

| CYP      | HLA    | Mitochondrial | Others |
|----------|--------|---------------|--------|
| CYP2D6   | HLA-A  | MTX1          | ABO1   |
| CYP2C19  | HLA-B  | MTX2          | ABO2   |
| CYP2C9   | HLA-C  | MTX3          | ABO3   |
| CYP2C10  | HLA-D  | MTX4          | ABO4   |
| CYP2C11  | HLA-E  | MTX5          | ABO5   |
| CYP2C12  | HLA-F  | MTX6          | ABO6   |
| CYP2C13  | HLA-G  | MTX7          | ABO7   |
| CYP2C14  | HLA-H  | MTX8          | ABO8   |
| CYP2C15  | HLA-I  | MTX9          | ABO9   |
| CYP2C16  | HLA-J  | MTX10         | ABO10  |
| CYP2C17  | HLA-K  | MTX11         | ABO11  |
| CYP2C18  | HLA-L  | MTX12         | ABO12  |
| CYP2C19  | HLA-M  | MTX13         | ABO13  |
| CYP2C20  | HLA-N  | MTX14         | ABO14  |
| CYP2C21  | HLA-O  | MTX15         | ABO15  |
| CYP2C22  | HLA-P  | MTX16         | ABO16  |
| CYP2C23  | HLA-Q  | MTX17         | ABO17  |
| CYP2C24  | HLA-R  | MTX18         | ABO18  |
| CYP2C25  | HLA-S  | MTX19         | ABO19  |
| CYP2C26  | HLA-T  | MTX20         | ABO20  |
| CYP2C27  | HLA-U  | MTX21         | ABO21  |
| CYP2C28  | HLA-V  | MTX22         | ABO22  |
| CYP2C29  | HLA-W  | MTX23         | ABO23  |
| CYP2C30  | HLA-X  | MTX24         | ABO24  |
| CYP2C31  | HLA-Y  | MTX25         | ABO25  |
| CYP2C32  | HLA-Z  | MTX26         | ABO26  |
| CYP2C33  | HLA-AA | MTX27         | ABO27  |
| CYP2C34  | HLA-AB | MTX28         | ABO28  |
| CYP2C35  | HLA-AC | MTX29         | ABO29  |
| CYP2C36  | HLA-AD | MTX30         | ABO30  |
| CYP2C37  | HLA-AE | MTX31         | ABO31  |
| CYP2C38  | HLA-AF | MTX32         | ABO32  |
| CYP2C39  | HLA-AG | MTX33         | ABO33  |
| CYP2C40  | HLA-AH | MTX34         | ABO34  |
| CYP2C41  | HLA-AI | MTX35         | ABO35  |
| CYP2C42  | HLA-AJ | MTX36         | ABO36  |
| CYP2C43  | HLA-AK | MTX37         | ABO37  |
| CYP2C44  | HLA-AL | MTX38         | ABO38  |
| CYP2C45  | HLA-AM | MTX39         | ABO39  |
| CYP2C46  | HLA-AN | MTX40         | ABO40  |
| CYP2C47  | HLA-AO | MTX41         | ABO41  |
| CYP2C48  | HLA-AP | MTX42         | ABO42  |
| CYP2C49  | HLA-AQ | MTX43         | ABO43  |
| CYP2C50  | HLA-AR | MTX44         | ABO44  |
| CYP2C51  | HLA-AS | MTX45         | ABO45  |
| CYP2C52  | HLA-AT | MTX46         | ABO46  |
| CYP2C53  | HLA-AU | MTX47         | ABO47  |
| CYP2C54  | HLA-AV | MTX48         | ABO48  |
| CYP2C55  | HLA-AW | MTX49         | ABO49  |
| CYP2C56  | HLA-AX | MTX50         | ABO50  |
| CYP2C57  | HLA-AY | MTX51         | ABO51  |
| CYP2C58  | HLA-AZ | MTX52         | ABO52  |
| CYP2C59  | HLA-AA | MTX53         | ABO53  |
| CYP2C60  | HLA-AB | MTX54         | ABO54  |
| CYP2C61  | HLA-AC | MTX55         | ABO55  |
| CYP2C62  | HLA-AD | MTX56         | ABO56  |
| CYP2C63  | HLA-AE | MTX57         | ABO57  |
| CYP2C64  | HLA-AF | MTX58         | ABO58  |
| CYP2C65  | HLA-AG | MTX59         | ABO59  |
| CYP2C66  | HLA-AH | MTX60         | ABO60  |
| CYP2C67  | HLA-AI | MTX61         | ABO61  |
| CYP2C68  | HLA-AJ | MTX62         | ABO62  |
| CYP2C69  | HLA-AK | MTX63         | ABO63  |
| CYP2C70  | HLA-AL | MTX64         | ABO64  |
| CYP2C71  | HLA-AM | MTX65         | ABO65  |
| CYP2C72  | HLA-AN | MTX66         | ABO66  |
| CYP2C73  | HLA-AO | MTX67         | ABO67  |
| CYP2C74  | HLA-AP | MTX68         | ABO68  |
| CYP2C75  | HLA-AQ | MTX69         | ABO69  |
| CYP2C76  | HLA-AR | MTX70         | ABO70  |
| CYP2C77  | HLA-AS | MTX71         | ABO71  |
| CYP2C78  | HLA-AT | MTX72         | ABO72  |
| CYP2C79  | HLA-AU | MTX73         | ABO73  |
| CYP2C80  | HLA-AV | MTX74         | ABO74  |
| CYP2C81  | HLA-AW | MTX75         | ABO75  |
| CYP2C82  | HLA-AX | MTX76         | ABO76  |
| CYP2C83  | HLA-AY | MTX77         | ABO77  |
| CYP2C84  | HLA-AZ | MTX78         | ABO78  |
| CYP2C85  | HLA-AA | MTX79         | ABO79  |
| CYP2C86  | HLA-AB | MTX80         | ABO80  |
| CYP2C87  | HLA-AC | MTX81         | ABO81  |
| CYP2C88  | HLA-AD | MTX82         | ABO82  |
| CYP2C89  | HLA-AE | MTX83         | ABO83  |
| CYP2C90  | HLA-AF | MTX84         | ABO84  |
| CYP2C91  | HLA-AG | MTX85         | ABO85  |
| CYP2C92  | HLA-AH | MTX86         | ABO86  |
| CYP2C93  | HLA-AI | MTX87         | ABO87  |
| CYP2C94  | HLA-AJ | MTX88         | ABO88  |
| CYP2C95  | HLA-AK | MTX89         | ABO89  |
| CYP2C96  | HLA-AL | MTX90         | ABO90  |
| CYP2C97  | HLA-AM | MTX91         | ABO91  |
| CYP2C98  | HLA-AN | MTX92         | ABO92  |
| CYP2C99  | HLA-AO | MTX93         | ABO93  |
| CYP2C100 | HLA-AP | MTX94         | ABO94  |
| CYP2C101 | HLA-AQ | MTX95         | ABO95  |
| CYP2C102 | HLA-AR | MTX96         | ABO96  |
| CYP2C103 | HLA-AS | MTX97         | ABO97  |
| CYP2C104 | HLA-AT | MTX98         | ABO98  |
| CYP2C105 | HLA-AU | MTX99         | ABO99  |
| CYP2C106 | HLA-AV | MTX100        | ABO100 |

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